TECHNICAL INFORMATION

**2024/08 (ver. 6)

*2024/07 (ver. 5)

Approval number: 30200BZI00022000

Category: Program 01 Disease diagnostic program

Highly Controlled Medical Device

Somatic Cell Mutation Analysis Program (for anti-cancer drug adaptation judgment) (70159013)

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【Shape, Structure and Principle, etc.】

Description for diagnostic tests using this system (device)

MyChoice診断システム is a companion CDx program which is a next generation sequencing-based *in vitro* diagnostic test, that assesses homologous recombination deficiency through tumor genomic instability (GI) [via an algorithmic measurement of Loss of Heterozygosity (LOH), Telomeric Allelic Imbalance (TAI), and Large-scale State Transitions (LST)] and detects and classifies sequence variants and large rearrangements (LR) in the *BRCA1* and *BRCA2* genes, using DNA obtained from fixed tumor tissue.

The patient specimen is packaged together with a test request form and held for pickup by a designated distributor in Japan. The specimen is collected and shipped by the distributor to the Myriad laboratory in Salt Lake City, UT, USA.

Each variant identified is checked against a database of previously classified variants. Variants not previously classified in the variant database, are classified by applying a defined set of objective criteria from multiple sources. The process and criteria (evidence) used to classify variants have been established in accordance with the recommendations of the American College of Medical Genetics and Genomics (ACMG) for standards in the interpretation and reporting of sequence variations (Richards et al 2015) and are consistently applied to all variants identified by Myriad. *BRCA1* and *BRCA2* variants are classified into one of the following five categories:

- Deleterious
- Suspected Deleterious
- Variant of Uncertain Significance
- Favor Polymorphism
- Polymorphism

Clinically significant variants are those classified as “deleterious” or “suspected deleterious”. Myriad’s variant database contains over 19,000 classified *BRCA1* and *BRCA2* variants. Less than 1% of the variants that will be identified during the sequencing process will be novel, and therefore require review for classification. Variants seen for the first time will be sent to the New Mutations Committee, which is composed of Myriad employees, for classification.

This product’s diagnostic program provides a Health Care Provider (HCP) with the combined interpretation of the sequencing and large rearrangement (LR) analyses and the Genomic Instability

Status (GIS). A positive Myriad HRD result due to either the presence of a pathogenic mutation in *BRCA1* or *BRCA2* (sequencing or LR) or a GIS above a defined threshold is correlated with homologous recombination deficiency (HRD).

【Intended Use or Indication】

This device is used to detect homologous recombination deficiency (HRD) by assessing Genomic Instability Status (GIS) and *BRCA1* and *BRCA2* gene mutations in genomic DNA extracted from tumor specimens. Results are used as an aid to determine the eligibility of patients with ovarian cancer for niraparib tosylate hydrate monotherapy treatment or olaparib in combination with bevacizumab treatment. Eligibility of patients for olaparib monotherapy treatment should be based on *BRCA1* and *BRCA2* gene mutations only.

****【Usage Method, etc.】**

[Instructions for Use with MailGate]

After testing has been completed, a secure email portal will be used for viewing the list of variants identified through testing, the Genomic Instability score, and the patient's Test Report. An email notification containing a link to a secure email portal will be sent to the Health Care Provider (HCP) indicated on the test request form. Within the secure email portal will be a secure email from Myriad that will provide the Interim Report, which contains a list of all variants identified from testing and the Genomic Instability score.

The HCP will be required to follow the following steps:

- 1) Click on the link to the secure email portal provided in the email
- 2) Enter the required login information (Registration for access to the secure email portal will be required when the secure email portal is accessed for the first time)
- 3) Open the secure email from Myriad to open the Interim Report
- 4) Review the list of variants and Genomic Instability score and request the patient's Test Report by selecting, "View Final Test Report"
- 5) View and print the patient's Test Report

[Instructions for Use with the Myriad Provider Portal]

After testing has been completed, a secure portal will be used for viewing the patient's Interim Report, a list of variants identified through testing and the Genomic Instability score, as well as the patient's Test Report. An email notification containing a link to a secure portal will be sent to the Health Care Provider (HCP) indicated on the test request form. Within the secure portal, the HCP can access the patient's Interim Report and Test Report.

The HCP will be required to follow the following steps:

- 1) Click on the link to the secure portal provided in the email.
- 2) Enter the required login information (registration will be required prior to initial login to access the secure portal).
- 3) Review the patient's Interim Report.
- 4) Click on the “承認” button to accept the Interim Report.
- 5) Click “結果を表示” to request the patient's Test Report.

6) View and print the patient's Test Report.

The Interim Report contains the list of all variants identified through the testing of the patient's sample, except for those variants for which definitive evidence indicates a lack of clinical significance (i.e., Favor Polymorphisms and Polymorphisms). The patient's Test Report will include the Genomic Instability Status (GIS) and the Tumor Mutation *BRCA1/BRCA2* Status for the patient. A positive or negative GIS is determined by a numeric threshold, based on the patient's Genomic Instability score (an algorithmic measurement of LOH, TAI, and LST). The classification and information for each variant that is clinically significant (deleterious or suspected deleterious) will be listed on the patient report. When multiple variants are provided on the report, the overall interpretation of the Tumor Mutation *BRCA1/BRCA2* Status is based on the variant with the highest clinical significance. Variants of uncertain significance and other findings where the clinical significance is unknown are listed on the second page of the patient report.

Patients evaluated with the MyChoice診断システム test that are determined to have a positive Genomic Instability Status or to carry a deleterious or suspected deleterious *BRCA1* or *BRCA2* mutation can be considered for treatment with niraparib tosilate hydrate or olaparib in combination with bevacizumab under the supervision of a HCP. Patients evaluated with the MyChoice診断システム test that are determined to carry a deleterious or suspected deleterious *BRCA1* or *BRCA2* mutation can be considered for treatment with olaparib monotherapy under the supervision of a HCP. It is recommended that these test results be communicated to the patient in a setting that includes appropriate counselling due to potential implications for the patient and family members. Variants of Uncertain Significance (VUS) are variants whose clinical significance has not yet been determined. Patients determined to carry a VUS cannot be considered for treatment with niraparib tosilate hydrate or olaparib if they do not also carry a deleterious or suspected deleterious *BRCA1* or *BRCA2* mutation.

A new sample may be requested if deemed appropriate. A Customer Service representative will contact the HCP to request an additional specimen and explain the reason a new specimen is being requested.

【Precautions for Use】

<Important Precautions>

- It should be confirmed that medical sites satisfy the conditions specified by the related societies before requesting a genetic test.
- Eligibility for treatment with drug should be determined after the patient report has been reviewed in its entirety. More specifically, a decision for eligibility for treatment with drug should not be made based solely on the results of the "Myriad HRD Status" at the beginning of the patient report, but after consideration of both the Genomic Instability Status and/or Tumor Mutation *BRCA1/BRCA2* Status. The decision to prescribe niraparib tosilate hydrate or olaparib in combination with bevacizumab should be based on the Genomic Instability Status and/or Tumor Mutation *BRCA1/BRCA2* Status (i.e., positive Myriad HRD Status). The decision to prescribe olaparib monotherapy should be based only on the Tumor Mutation *BRCA1/BRCA2* Status.
- The test is designed to detect tumor DNA variants within specific regions of the *BRCA1* and *BRCA2* genes. This test does not determine whether variants are germline or somatic. Gene coding regions and portions of non-coding intronic regions are analyzed by sequence analysis and typically do not extend more than 20 base pairs (bp) proximal to the 5' end and 10 bp distal to the 3' end of each exon. Additional variants outside of the assessed regions that may impact patient care will not be detected.

- Unequal allele enrichment may result from long indels or rearrangements under probe. Indels > 25 bp in length can be detected by this assay. However, the ability to detect any particular indel may be impacted by the location and nature of the mutation, the local sequence context, the DNA quality, and the non-tumor DNA content in the sample provided.
- The test has been designed to detect genomic rearrangements (i.e., deletions or duplications) involving the promoter and coding exons of *BRCA1* and *BRCA2*, however, the detection of large rearrangement deletions and duplications is dependent on the quality of the submitted specimen. Whole gene duplications and deletions may not be detected by the MyChoice 診断システム assay. Other terminal duplications are reported as variants of uncertain significance.
- There may be uncommon genetic abnormalities such as specific insertions, inversions, and certain regulatory mutations that will not be detected by this product. This analysis, however, is believed to rule out the majority of abnormalities in the genes analyzed.
- The classification and interpretation of all variants identified reflects the current state of scientific understanding at the time the result report is issued. In some instances, the classification and interpretation of variants may change as scientific information becomes available.
- This system is designed to detect sequencing variants and genomic rearrangements (deletions or duplications) involving the promoter and coding exons of *BRCA1* and *BRCA2*. As such, use in the following cases may lead to false negative results.
 - When there are variants that may impact patient care outside of the assessed regions.
 - When mutations result in some types of errors in RNA transcript processing
 - When mutations are insertions that do not result in duplications
- Alterations at allele frequencies below the established limit of detection may not be detected consistently. The sequence variant allele frequency LOD is ~ 6.5-8% for samples run at 200ng and ~10-20% for samples run at 7.5ng. The LR allele frequency LOD is ~20-30% for samples run at 200ng and ~40-50% for samples run at 7.5ng.
- The Genomic Instability Status (GIS) may not be determined if the ratio of tumor to benign nuclei of the specimen is below the established LOD. At optimal DNA input levels (e.g., 200ng) the GIS LOD is ~ 20%-30% tumor cell percentage. At the lowest DNA input level (7.5ng) the LOD is ~40-49% tumor cell percentage.

【Other Precautions】

- When this product is used, the latest Electronic Technical Information for niraparib tosilate hydrate or olaparib in Japan should be referenced.
- This assay may provide information on the risk of Hereditary Breast and Ovarian Cancer syndrome. However, this assay has not been approved for such a purpose, and the Japanese regulatory authority has not acknowledged that its clinical significance or analytic validity are appropriate.

【Clinical Significance】

I. Summary of Clinical Study - Niraparib Tosilate Hydrate PR-30-5020-C (QUADRA)

The niraparib tosilate hydrate clinical study PR-30-5020-C (QUADRA) was a multicenter, open-label, single-arm clinical trial designed to evaluate the safety and efficacy of niraparib tosilate

hydrate in patients with advanced, relapsed, high-grade serous epithelial ovarian, fallopian tube or primary peritoneal cancer who had received three or more previous chemotherapy regimens.

a) Accountability of Premarket Approval (PMA)

To determine patient eligibility, a tumor sample was sent to a centralized laboratory for immediate myChoice HRD CDx (HRD) testing. Archival or fresh tumor tissue was required. *(Note: for patients enrolled after implementation of protocol version 3 [amendment 2; 24 May 2016], the sample was required prior to enrollment and could be sent in advance of the protocol-defined screening period in order to facilitate the screening and enrollment process).* Patients were to wait for the results from the on-study centralized HRD testing prior to enrollment unless they had a previously detected germline breast cancer gene (*BRCA*) mutation (*gBRCAm*). Blood samples were also collected for all patients during screening for determination of *gBRCAm* status. If *gBRCA* status was previously confirmed positive, then it was not necessary to wait for HRD testing results for enrollment into the study; however, confirmatory HRD testing still needed to be performed.

b) Effectiveness Results

The primary efficacy endpoint was Investigator-assessed confirmed objective response rate (ORR) as defined by RECIST (version 1.1). The primary analysis population for this endpoint included all patients with HRD positive tumors who received 3 or 4 prior lines of therapy (LOT) and whose disease was platinum-sensitive to the last platinum based therapy. Patients with prior poly (ADP-ribose) polymerase (PARP) inhibitor (PARPi) treatment were not included in the primary analysis population. The patient's best overall response (BOR) was determined based on the overall responses at all timepoints between the date of the first dose and the date of first documented radiological disease progression, the date of subsequent anticancer therapy, or the date of study discontinuation, whichever occurred first. Patients with a BOR of either confirmed complete response (CR) or confirmed partial response (PR) were considered to have responded to treatment ("responders"). All other patients were considered not to have responded to treatment ("non-responders"). Post hoc analyses were also conducted in subgroups defined by platinum sensitivity and biomarkers.

Investigator-assessed confirmed ORR as defined by RECIST (version 1.1) was analyzed for the biomarker-defined population overall and in the following biomarker subsets:

- *tBRCAm* (regardless of platinum sensitivity)
- Non-*tBRCA*/GIS positive platinum-sensitive

In the biomarker-defined population, meaningful ORR and duration of response (DOR) were observed among all subgroups, including the following:

- In the HRD positive cohort (n=98), the ORR was 24.4% and the median DOR was 8.3 months.
- In patients with *tBRCAm* tumors (n=63), the ORR was 28.6% and median DOR was 9.2 months. When this group was analyzed by platinum sensitivity status, patients with *tBRCAm* platinum-sensitive disease (n=18) had an ORR of 38.9%; patients with *tBRCAm* platinum-resistant disease (n=21) had an ORR of 28.5%; and patients with *tBRCAm* platinum-refractory disease (n=16) had an ORR of 18.8%.

- In patients with non-*tBRCA*/GIS positive platinum-sensitive disease (n=35), the ORR was 20.0% and median DOR was 6.6 months.

c) Conclusions

In summary, the efficacy results from the PR-30-5020-C (QUADRA) clinical study demonstrated that the magnitude of clinical benefit expected from niraparib tosilate hydrate treatment can be predicted by the clinical benefit continuum defined by clinical and molecular biomarkers. In line with the continuum of benefit based on molecular and clinical biomarkers, the highest response rates were observed among patients with tumors that were GIS positive and platinum-sensitive or *tBRCA*m regardless of platinum sensitivity status.

**II. Summary of Clinical Study - Olaparib D0817C00003 (PAOLA-1)

The clinical utility of the Myriad myChoice[®] CDx has been evaluated using ovarian cancer patient samples collected from the PAOLA-1 clinical trial.

1. Overview of PAOLA-1 Study

A phase III, randomized, double-blind, placebo-controlled, multicenter study was conducted in 806 patients (537 patients in the olaparib/bevacizumab group, 269 patients with placebo/bevacizumab) with newly diagnosed, advanced (FIGO stage IIIB-IV), high grade epithelial ovarian, fallopian tube or primary peritoneal cancer treated with standard first-line treatment, combining platinum-taxane chemotherapy and bevacizumab concurrent with chemotherapy and in maintenance to evaluate the efficacy and safety of olaparib (tablet) 300 mg given twice daily in combination with bevacizumab) in comparison with placebo in combination with bevacizumab.

Additionally, PFS was investigated in exploratory biomarkers subgroups using the Myriad myChoice[®] Plus test to further understand the consistency of treatment effect of olaparib across potential predictive and prognostic factors.

Myriad HRD status is based on the Genomic Instability Status (GIS) and/or Tumor Mutation *BRCA1/BRCA2* (*tBRCA*m) Status. A positive Myriad HRD status is determined either by presence of a *tBRCA1/2* mutation or by a Genomic Instability (GI) score at or above a cut-off of 42 in the absence of a *BRCA1/2* mutation.

Of the 806 randomised patients in PAOLA-1, 755 (93.7%) had a Myriad *tBRCA* available and 664 (82.4%) had an available Myriad HRD status with the myChoice[®] HRD plus RUO test. With the myChoice[®] CDx, a total of 755/806 (93.7%) patients were considered for testing. A total of 700/806 (86.8%) patients had a valid *tBRCA* result and 643/806 (79.8%) had a valid HRD status result.

2. Concordance between the myChoice[®] HRD plus RUO test and the myChoice[®] CDx in PAOLA-1

The concordance between the Myriad myChoice[®] HRD plus RUO test and Myriad's myChoice[®] CDx for HRD status at cut-off 42 for patients with a valid result for both tests, was as follows:

OPA: 97.8% (95% CI 96.3%, 98.8%)

PPA: 98.6% (95% CI 96.8%, 99.6%)

NPA: 96.6% (95% CI 93.7%, 98.4%)

3. Comparison of efficacy between PAOLA-1 full analysis set and the subgroups determined by Myriad myChoice® CDx

PAOLA-1 met its primary objective, demonstrating a statistically significant and clinically meaningful improvement in PFS in the FAS for olaparib vs placebo when added to bevacizumab (HR 0.59; 95% CI 0.49 to 0.72; p<0.0001). Olaparib prolonged PFS by a median of 5.5 months over placebo. The clinical outcome data for the Full Analysis Set and the efficacy data from post-database lock exploratory analyses performed in the subset of patients for whom Myriad HRD status for cut-off ≥42 defined by the Myriad myChoice® CDx test are shown in the following table.

Clinical outcome (PFS by investigator assessment) of PAOLA-1

FAS		myChoice® CDx Myriad HRD positive		myChoice® CDx Myriad HRD status negative		myChoice® CDx Myriad HRD status unknown	
Olaparib/bevacizumab (n=537)	Placebo/bevacizumab (n=269)	Olaparib/bevacizumab (n=251)	Placebo/bevacizumab (n=125)	Olaparib/bevacizumab (n=180)	Placebo/bevacizumab (n=87)	Olaparib/bevacizumab (n=106)	Placebo/bevacizumab (n=57)
Number of events/total number of patients (%)							
280/537 (52.1)	194/269 (72.1)	87/251 (34.7)	86/125 (68.8)	135/180 (75.0)	68/87 (78.2)	58/106 (54.7)	40/57 (70.2)
Median PFS (months)							
22.1	16.6	39.3	17.7	16.6	16.2	22.1	15.3
HR (95% CI)							
0.59 (0.49-0.72)		0.35 (0.26-0.48)		1.00 (0.75-1.34)		0.68 (0.45-1.02)	

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III. Summary of Clinical Study - Olaparib D0818C00001 (SOLO1)

The clinical utility of the Myriad myChoice® CDx has been evaluated using ovarian cancer patient samples collected from the SOLO1 clinical trial.

1. Overview of SOLO1 Study

SOLO1 (Study Code D0818C00001) is a Phase III, randomised, double-blind, placebo-controlled, multi-centre study to assess the efficacy of olaparib maintenance monotherapy in high risk advanced ovarian cancer patients (including patients with primary peritoneal and / or fallopian tube cancer) with a documented mutation in *BRCA1* or *BRCA2* that is predicted to be deleterious or suspected deleterious (known or predicted to be detrimental/lead to loss of function) who have responded following first line platinum based chemotherapy.

The SOLO1 trial was designed to recruit *BRCAm* patients i.e., germline *BRCAm* (*gBRCAm*) or somatic *BRCAm* (*sBRCAm*). The entry criteria for patients in SOLO1 included the requirement to have a deleterious or suspected deleterious *BRCA1* or *BRCA2* mutation (from blood or tumour specimens), determined prior to study entry from an existing local BRCA mutation test or from prospective testing using the Myriad BRCAAnalysis® or BRCAAnalysis CDx test. Local results

that could not be confirmed by the Myriad assay were excluded from the final analysis. Due to sample export restrictions, local and prospective testing for patients within China was performed by a laboratory within the country.

SOLO1 met its primary objective, demonstrating a substantial improvement in progression-free survival (PFS) by investigator assessment for patients treated with olaparib compared with placebo. This PFS improvement was statistically significant and clinically relevant, as evidenced by the magnitude of effect: a 70% reduction in the risk of disease progression or death at any point in time for olaparib vs placebo-treated patients (HR 0.30; 95%CI 0.23 to 0.41; $p < 0.0001$), with a median PFS not reached on the olaparib arm vs 13.8 months for placebo.

2. Concordance between Myriad *tBRCA* test and Myriad *gBRCA* test in SOLO1

It is to be expected that a tumour *BRCA* test and a germline *BRCA* test will not be fully concordant since germline tests are unable to detect somatic variants. In addition, there may be uncommon genetic abnormalities such as specific insertions, inversions, and certain regulatory mutations that will not be detected. Overall, treatment eligibility was concordant in 284/292 (97.2%, 8 patients reported as non-*tBRCAm*)

- Percent positive agreement between the Myriad *gBRCA* and Myriad *tBRCA* tests results was calculated as 97.6% (95% CI 95.07% to 99.02%).
- Negative percent agreement between cases confirmed as germline *BRCA* negative and a Myriad *tBRCA* test have been presented for LIGHT study. The negative percent agreement was 88.8% (95% CI, 84.3%, 92.3%).

3. Comparison of efficacy between SOLO1 full analysis set, Myriad *tBRCAm* subset, Myriad *gBRCAm* subset

Post-database lock exploratory analyses indicated that efficacy results are consistent between the FAS, the Myriad *tBRCAm* subset defined by the myChoice[®] CDx test and Myriad *gBRCAm* subset defined by the BRACAnalysis[®] CDx test.

Table Summary of key efficacy outcome variables (Full Analysis Set, Myriad *tBRCAm* subset, Myriad *gBRCAm* subset)

	Full Analysis Set n=391		Myriad <i>tBRCA</i> n=284		Myriad <i>gBRCA</i> n=383	
	Olaparib 300 mg bd (n=260)	Placebo (n=131)	Olaparib 300 mg bd (n=191)	Placebo (n=93)	Olaparib 300 mg bd (n=253)	Placebo (n=130)
PFS by Investigator Assessment (maturity 50.6%)						
Number of events/total number of patients (%)	102/260 (39.2)	96/131 (73.3)	75/191 (39.3)	69/93 (74.2)	99/253 (39.1)	95/130 (73.1)
Median PFS (months)	Not reached	13.8	Not reached	13.8	Not reached	13.8
HR (95% CI)	0.30 (0.23-0.41)		0.29 (0.21-0.41)		0.30 (0.22-0.40)	
P-value (2-sided)	p<0.0001		p<0.0001		p<0.0001	

【Approval Conditions】

1. It is required to undergo the document-based or on-site inspection in compliance with the provisions of Article 23-2-5, Paragraph 9 of the Pharmaceutical and Medical Device Act (Act No.145 of 1960, hereinafter referred to as the “Act”) pursuant to the provisions of Article 23-2-17, Paragraph 5 of the Act, and to submit an operational update report of the variant classification process to PMDA one year after approval and annually thereafter for the time being.
2. The applicant is required to take necessary measures and appropriate management strategies to ensure that tumor samples and information obtained from them are not used for the purposes other than the items specified in the attached application form, and also implement the latest measures to ensure security and privacy to prevent unauthorized access.
3. The quality control for information of variant classification and genomic instability must be conducted as specified in the Remarks section of the attached application. When it is planned to revise the quality control method for information of variant classification and genomic instability specified in the Remarks section of the attached application (except for a minor revision defined in the MHLW Ordinance specified in Article 23-2-5, Paragraph 15 of the Act), the revision must be approved by the Minister of Health Labor and Welfare in compliance with the provisions of Article 23-2-5, Paragraph 15 of the Act pursuant to the provisions of Article 23-2-17, Paragraph 5 of the Act. The applicant should note that, for the approval of the above revision, the provisions of Article 23-2-5, Paragraph 17, Article 23-2-6 and Article 23-2-7 of the Act are applied pursuant to the provisions of Article 23-2-17, Paragraph 6 of the Act.

***【Name or Designation etc. of Marketing Authorization Holder or Manufacturer】**

Foreign Special Approval Holder

Myriad Genetic Laboratories, Inc. (*in Japanese*)

(Myriad Genetic Laboratories, Inc.)

US

DMAH

Myriad Genetics G.K. (*in Japanese*)

(Myriad Genetics G.K.)

TEL: 03-4520-7080

Foreign Manufacturer

Myriad Genetic Laboratories, Inc. (*in Japanese*)

(Myriad Genetic Laboratories, Inc.)

US

Myriad Genetic Laboratories, Inc., MyChoice診断システム[®] and their respective logos, are either trademarks or registered trademarks of Myriad Genetics, Inc. and its subsidiaries in the United States and other jurisdictions.

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