

# Male or Female?

Get a SneakPeek.



**A >99% accurate early fetal sex DNA test, starting at six weeks into pregnancy<sup>1</sup>**

Simple blood test | Fast results

**Ask about SneakPeek Clinical today!**

1. Szucs HD, Asprer L, Abunadi N, Milot H, Cacia S, Jacob C (2022) SneakPeek Early Gender Test: The Earliest and Most Accurate Method for Fetal Sex Determination at 6-Weeks Gestation. Clinics Mother Child Health. 19:419.

## How does SneakPeek® work?

While every person has their own DNA in their bloodstream, a pregnant person's also contains DNA from their unborn child. SneakPeek looks for male chromosomes in a small sample of a patient's blood. If male chromosomes are found, the baby is predicted to be male. If none are found, the baby is predicted to be female. SneakPeek is >99% accurate at six weeks into pregnancy.<sup>1</sup>



## Learn more about SneakPeek Clinical

Scan the QR code



Myriad Genetics / 322 North 2200 West, Salt Lake City, UT 84116

©2023 Myriad Genetics, Inc. All rights reserved. SneakPeek is a registered trademark of Gateway Genomics, LLC. Myriad Genetics, Health. Illuminated., and its respective logos are either trademarks or registered trademarks of Myriad Genetics, Inc. and its subsidiaries in the United States and other jurisdictions. SPPATSEX 0723